




Nume Peptides Inc
info@numepeptides.com
http://www.numepeptides.com

SAMPLE NAME:
SAMPLE CODE:
RECEIVED DATE:
PUBLISHED DATE:

Melanotan II
7RQY-43PD
09/23/2026
09/25/2026



DECLARED IDENTITY: Melanotan II
MATRIX TYPE: PEPTIDE
SAMPLE NAME: Melanotan II
SAMPLE SIZE: 10.0 mg
LOT CODE: 070326




9.816 mg
Total Quantity



99.862%
Total Blend Purity

Sample	Test	Specification	Result	Status
Melanotan II	IDENTITY	Melanotan II	Melanotan II	CONFORMS
Melanotan II	QUANTITY	As measured	9.816 mg	MEASURED
Melanotan II	PURITY	≥ 98 %	99.862%	CONFORMS

HPLC Chromatogram Report



RESULTS REVIEWED BY:

Anthony Burke PhD
Lab Manager
09/25/2026

RESULTS CERTIFIED BY:

Dr. Levi Fried MD
Lab Director
09/25/2026

ACCUMARKLABS.COM
7RQY-43PD
09/25/2026





Nume Peptides Inc
info@numepeptides.com
http://www.numepeptides.com

SAMPLE NAME: Melanotan II
SAMPLE CODE: 7RQY-43PD
RECEIVED DATE: 09/23/2026
PUBLISHED DATE: 09/25/2026

Endotoxin USP85 LAL

Test	Specification (EU/mL)	Result (EU/mL)	LOQ (EU/mL)	Verdict
Endotoxin USP85 LAL	≤5	< LOQ	2.5	Conforms

Method: Endotoxin testing performed using Limulus Amebocyte Lysate assay in accordance with USP <85> under validated laboratory conditions.

Rapid Sterility Screening (PCR)

Test	Specification	Result	Unit	Verdict
Rapid Sterility Screening (PCR)	Not Detected	Not Detected		Conforms

Method: Rapid sterility screening performed by real-time polymerase chain reaction (qPCR) targeting conserved microbial DNA sequences under validated laboratory conditions.

Fentanyl Screening

Basis: Qualitative lateral-flow immunoassay | Strip detection threshold 200 ng/mL in prepared solution

Test	Specification	Result	Unit	Verdict
Fentanyl	Not Detected	Not Detected		Conforms

Method: Lateral-flow immunoassay test strip | Sample Prep: Aqueous reconstitution prepared ; test strip immersed directly per manufacturer procedure

Reporting Fentanyl screening is qualitative. Results are reported as Not Detected or Detected relative to the test strip's stated detection threshold; no numeric concentration is measured or implied. Not Detected indicates fentanyl was not present above that threshold in the prepared sample solution.

RESULTS CERTIFIED BY:


Anthony Burke PhD
Lab Manager
09/25/2026

RESULTS CERTIFIED BY:


Dr. Levi Fried MD
Lab Director
09/25/2026

ACCUMARKLABS.COM
7RQY-43PD
09/25/2026

